## PREMACRYL® PLUS

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en

## Specification:

PREMACRYL® PLUS is a modified self-curing methacrylate resin designed for orthodontic use. One of its features is the reduced dissolution time of the polymer beads which means that the dough maintains its plasticity for a longer period. An extended range of colour shades available for the resin is another new feature. The resin is now available in six shades: R – red, GR – green, G – yellow, B – blue, O – clear and V – pink. Unlike the traditional base resins, only the liquid is coloured with Premacryl® Plus, except for the O and V shades.

**Indications:**

- Manufacture of orthodontic apparatuses.
- Manufacture of immediate dentures.
- Repair of removable dentures.

**Contra-indications:**

The application of PREMACRYL® PLUS is contra-indicated in the case of known hypersensitivity of the patient to any of the components of the product. Undesirable side-effects: Irritation in the prosthetic bed area can be seen among patients sensitive to methylmethacrylate.

**Working:**

The mixture is prepared in a glass or ceramic jar. After mixing, the mixed resin must be left still, to avoid porosity which would be caused by further mixing. In the traditional procedure, the dough is adapted to the insulated plaster model. At this stage, the retention elements are anchored in the dough. The prepared orthodontic apparatus is polymerised in a simple pressure polymerisation device working on the principle of water column compression.

**Mixing ratio:**

*By weight* - 2 gr. of powder to 1 gr. of liquid

*By volume* - 3:1 (3 parts of powder to 1 part of liquid)

**Mixing time:** 1.5 min

**Working time:** 5–6 min. at a temperature of 25 °C (from the start of mixing to the start of setting)

**Setting time:** 15 minutes at a pressure of 0.2 – 0.3 MPa, in water at a temperature of 50–60 °C.

**Notice:**

It should be taken into consideration that the working procedure as well as the mixing are considerably accelerated if carried out at a higher temperature (23–25 °C).

**Warning:**

Liquid: H225 Highly flammable liquid and vapour. H315 Causes skin irritation. H335 May cause respiratory irritation. H317 May cause an allergic skin reaction.

P210 Keep away from heat/sparks/open flames/hot surfaces. – No smoking. P262 Do not get in eyes, on skin, or on clothing. P270 Do not eat, drink or smoke when using this product.

**Waste material disposal:**

The residue of unused product and contaminated packaging material should be transferred to a company authorized to dispose of waste materials. Observe your local legal regulations.

**Storage:** To be stored in a dry place at a temperature of 5–25 °C.

**Expiry period:** 4 years

**Package:**

*Pulvis:* 500 gr. – O clear and V pink

*Liquid:* 250 gr. – clear, red, blue, green and yellow

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The logo consists of the letters 'de' in a white, sans-serif font, centered within a black oval shape.**Charakteristik:**

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