## CZ

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Katétr Nelaton

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Foley Catheter Nelaton

#### DESCRIPTION

The Latex Foley Catheter is flexible tube that is passed through the urethra during urinary catheterization and into the bladder to drain urine, or for inserting fluids into the bladder. It consists of shaft, drainage funnel, inflation funnel, irrigation funnel (if present), balloon and valve. It is easy to use, reliable in performance, and irritation free.

The Latex Foley Catheter has the following features:

- Made from natural latex
- Good biocompatibility
- Silicone coated surface reducing allergic reaction
- Smooth tapered tip facilitating introducing
- Color-coded for visualization of size
- Alternative types are available
- sterile by EO
- For single use only

#### INTENDED USE/ INDICATIONS

Disposable urinary catheter is used to discharge the bladder, for short-term use. The inserted catheter is used to clamp urine by connecting it to a urinal bag.

#### [USE INFORMATION]

The Latex Foley Catheter is passed through the urethra during urinary catheterization and into the bladder to drain urine, or for inserting fluids into the bladder. It is used in departments of urology, internal medicine, surgery, obstetrics, and gynecology for drainage of urine and medication. It is also used for patients suffering from moving with difficulty or being completely bed-ridden.

#### CONTRAINDICATIONS

- Acute urethritis
- Acute prostatitis
- Acute epididymitis

#### WARNINGS / PRECAUTIONS

- Do not use if package has been opened or damaged.
- The leakage or deflation of balloon should be tested prior to use.
- The product is intended for use by physicians trained and aseptic technique must bepracticed.
- Do not over inflate the balloon. Refer to outer unit pack or funnel of catheter for ballooncapacity.
- nflate catheter balloon with sterile distilled water.
- Do not use ointment and lubricants having petroleum base. Use silicone oil or watersoluble lubricants.
- Do not clamp catheter shaft. It may damage catheter and prevent deflation.
- Do not use needle to inflate balloon. Use Luer syringe.
- Do not use after the expiry date.
- For single patient use only.
- This is a disposable product for single use, any resterilization or reuse of the product canresult in hazard to the patient.
- Precautions shall be taken when discarding device and disposal of the device shall bemad in accordance with applicable national regulations for biological hazardous waste.

#### ADVERSE REACTIONS

- Bleeding into or around the catheter
- Catheter draining little or no urine despite enough fluid intake
- Fever, chills
- Leakage of large amounts of urine around the catheter
- Urine with a strong smell or becomes thick or cloudy
- Urethral swelling around the catheter
- Urinary tract or kidney infections

#### •Urethritis

#### DIRECTIONS FOR USE

- Generously lubricate tip and shaft of the catheter. Do not use ointment and lubricantshaving petroleum base.
- Carefully insert catheter tip into the bladder (normally indicated by urine flow), and thena further 3cm to ensure balloon is also inside it.
- Using the syringe set inflate balloon. Recommended volume is marked on funnel ofcatheter.
- For deflation, cut off the inflation funnel above the valve, or using a syringe withoutneedle push into valve to facilitate drainage.
- The catheter cannot be retained in human body for over 14 days.

#### Caution:

- For single patient use only.
- Sterile if package is unopened or undamaged.
- The size and balloon capacity is marked on outer unit pack and funnel of catheter.
- This product contains latex and can trigger allergic reactions to people who have increased sensitivity to certain proteins found in latex.

#### STORAGE CONDITIONS

- Protect product from moisture and excessive heat.
- Avoid prolonged exposure to ultraviolet and fluorescent light.
- Store in manner preventing crushing.

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