CZ

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Instruction for use Spigot

Description

The spigot is used for hygienic stoppage of urine flow from catheters during nursing procedures. It is non-invasive device utilized for sealing the catheter for a short time and allow the urine to gather in the bladder.

The Spigot is made from PP. It has the following features:

- Made from non-toxic PP
- Universal size is formed to fit all sizes of catheters
- Graduated design makes it easy to accommodate the tubing of various internal diameters.
- Ribbed design provides strong and secure connections
- Ergonomic flanged design assists the withdrawal and maximal grip
- 100% latex-free
- Sterile pack-Sterilized by EO, single use

Intended use

The Spigot's intended use is to seal the drainage funnel of the Urethral Catheter as a prevention of nosocomial urinary tract infection.

Contraindications

- There are none known contraindications.

Warnings / precautions

- Do not use if package has been opened or damaged.
- Before use read all warnings and instruction for use, improper use can result in serious or fatal illness or injury.
- Do not use after expired expiry date.
- For single patient use only.
- This is a disposable product for single use, any re-sterilization or reuse of the product can result in hazard to the patient.
- Precautions shall be taken when discarding device and disposal of the device shall be met in accordance with applicable national regulations for biological hazardous waste.

Adverse reactions

- No adverse reaction has been reported.

Direction for use

- Clean the drainage funnel of the catheter using a sterile towel
- Take out one sterile packed spigot
- Open the sterile packing
- Insert the spigot into the drainage funnel of the catheter.

Caution:

- Single use only
- Sterile if packing is unopened or undamaged.
- Store in a dry, cool and dark place.

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CO., LTD.
Shanghai International Holding Corp. GmbH