<b>DE</b>	<b>Gebrauchsanweisung Dürr Universalkanülen / Chirurgische Kanülen</b>
<b>EN</b>	<b>Operating Instructions Dürr Universal Cannulas / Surgical cannula</b>
<b>FR</b>	<b>Instructions d'utilisation Canules universelles / Canules chirurgicales</b>
<b>ES</b>	<b>Instrucciones de uso Cánula Universales / Cánula Quirúrgicas</b>
<b>IT</b>	<b>Istruzioni d'uso Cannula universalì / Cannule chirurgiche</b>
<b>NL</b>	<b>Gebruiksaanwijzing Dürr universele canules /Chirurgische canule</b>
<b>PT</b>	<b>Instruções de uso Cânula universais / Cânula cirúrgicas</b>
<b>SV</b>	<b>Bruksanvisning Dürr Universalkanyler / Kirurgiska kanyler</b>
<b>DA</b>	<b>Brugsanvisning Dürr universalkanyle / kirurgisk sugekanyle</b>
<b>FI</b>	<b>Käyttöohjeet Dürr peruskanyylit / kirurgiset kanyylit</b>
<b>NO</b>	<b>Bruksanvisning Dürr Universalkanyler / kirurgisk kanyle</b>
<b>RU</b>	<b>Инструкция по эксплуатации и монтажу Дюр универсальные канюли / хирургические канюли</b>
<b>EL</b>	<b>Οδηγίες χρήσης Ρύγχη γενικής χρήσης Dürr / χειρουργικά ρύγχη</b>
<b>PL</b>	<b>Instrukcja użytkowania Uniwersalne kaniule Dürr / kaniula Chirurgiczne</b>
<b>CS</b>	<b>Návod k použití Dürr Univerzální kanyly / Chirurgické kanyly</b>
<b>SK</b>	<b>Návod na použitie Univerzálne kanyly Dürr / chirurgické kanyly</b>

9000-616-08

## Dürr Universalkanüle / Universalkanüle Petito / Chirurgische Kanüle

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## Dürr Universal Surgical / Universal Surgical Petito / Universal Surgical Cannula

A variety of cannula for the suction of secretions and spray mist from the mouth of the patient:

Universal Surgical - for adults

Universal Surgical Petito - for children

Universal Surgical Cannula\* - for suction during surgical treatment

### Preparation

Cannula are medical instruments which come into contact with mucous membranes and skin or tissue which has been affected by illness.

Remember to observe all current local and national regulations concerning the cleaning, disinfection, sterilisation and storage of cannula!

In the Federal Republic of Germany, for example, a risk evaluation must be carried out and medical products are to be categorised according to the Recommendations of the Commission on Hospital Hygiene and the Prevention of Infections. This medical product is categorised in group: "Semicritical/critical Group B with increased requirements for preparation and use".

### Transport / Storage

After use cannula must be placed in a suitable closed container which is filled with an instrument disinfectant solution (e.g. ID 212, ID 212 forte, ID 213; observe the manufacturer's instructions), collected and transported in a securely sealed vessel to the site of cleaning.

### Preparation for cleaning

No special requirements are necessary. Where necessary, remove larger dirt particles etc. using a disinfectant disposable wipe (e.g. FD 350).

### Cleaning and disinfection, automatic

Use cleaning and disinfection devices strictly according to manufacturer's instructions.

- Place the cannula into the unit in such a way that the interior surfaces will be rinsed thoroughly and the water can drain off.
- Set up the cleaning and disinfection cycles according to the manufacturer's instructions.
- After the cleaning and disinfection cycles have finished check the cannula carefully for signs of residual dirt and repeat the cycles if necessary.

### Cleaning and disinfection, manual

Cleaning takes place in a wet chemical process using a VAH/DGHM listed instrument disinfection agent:

- Place cannula fully submerged in an instrument-disinfectant and cleaning agent (e.g. Dürr ID 212, ID 212 forte; ID 213, observe manufacturer's instructions carefully)
- After the cleaning and disinfection cycles have finished check the cannula carefully for signs of residual dirt and repeat the cycles if necessary.
- Rinse cannula thoroughly with water.
- Drying should be carried out using hygienically clean, disposable wipes, or under the application of compressed air.

## Check and test of function

Frequent application and cleaning has little effect on cannula. The end of the product's life cycle is determined by the amount of wear and tear or damage during usage. Check the cannula visually and where necessary, damaged cannula must be disposed of. Before disposal cannula must be disinfected and sterilised.

## Packaging

Cannula can be sterilised either unpacked\* or packed in a container or in sterile product packaging. Use a suitable sterile product packaging!

## Steam Sterilisation

### WARNING

**The wrong sterilisation can reduce effectiveness and can lead to damage to the product**

- Only steam sterilisation is permissible.
- Do not exceed 134 °C.
- Observe the process parameters.
- Do not use any other process.
- Place filter paper on sterilisation tray before use.
- Always place objects for sterilisation together with instruments in order to achieve a balanced temperature distribution during sterilisation.
- Avoid metal contact between any of the objects for sterilisation and any instruments or with the tray edge itself.
- Place the objects to be sterilised in such a way that they do not come too close to the heating elements and do not come into contact with the internal walls and place on the middle tray.
- Only sterilise the product in a disassembled state.

Temperature: 134 °C

Gauge pressure above ambient environment: 2.16 bar / 0.216 MPa

Holding time: 5 min

- Sterilised implements must be sterilised in a small steam sterilisor before use according to European guidelines EN 13060 in sterilisation cycle B and subsequently dried.
- Observe the manufacturer's instructions concerning usage of the steam sterilisation unit and the correct placement of objects to be sterilised.
- Observe all national and local guidelines and regulations concerning sterilisation of medical products.

## Storage

Store the cannula so they are protected against any possible contamination.

## Important information

The instructions described above have been validated by Dürr Dental as suitable and are recommended for the preparation of cannula before and after use. The surgery manager is responsible for ensuring that the correct and safe preparation of the instruments is carried out using suitable appliances, materials and trained personnel in order to produce the desired results. Validation and routine monitoring of the preparation procedures are necessary. Any and every deviation from the procedures described above by preparation staff is solely the responsibility of the personnel in terms of effectivity and the possible negative repercussions.

## Canules universelles de Dürr / Canules universelles Petito / Canules chirurgicales

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DÜRR DENTAL AG